



**The powered
loupe for
unbound
performance.**



DRAGONFLY™
ORASCOPTIC

The era of the powered loupe has arrived.

Dragonfly™ by Orascoptic™

is a loupe with fully integrated battery power and circuitry to operate a built-in headlight.

A triumph of optical engineering, Dragonfly seamlessly combines Orascoptic's industry-leading optics and award-winning headlight into a single, lightweight product. All electronics and batteries are completely contained within the custom-built frame itself, so there are no visible wires or cables. An innovative counter-balance design affords all-day comfort by minimizing pressure on the user's nose bridge.

- Untethered freedom of movement.
- Magnification and illumination working in harmony.
- Balanced ergonomic comfort.
Reduce the weight on your nose by as much as 30% compared to traditional loupes.

Our vision is to drive the best patient outcomes by providing **Superior Visualization** for every healthcare procedure.

Challenge the status quo, elevate dental and medical procedures.



Our premium PRO model supports the full line of Orascoptic telescopes and is loaded with bonus features including extra batteries, multiple light intensity settings, and extremely durable ceramic paint finish.*



MAGNIFICATION

Traditional: 2.5x, 3.0x, 3.5x, 4.5x, 5.5x
Ergo: 3.5x, 5.0x
Variable: 2.5x-3.5x, 3x-4x-5x, 3.5x-6.5x


LIGHT INTENSITY

High: 78 Lumens
Low: 49 Lumens

BATTERY RUN-TIME (Set of 2)

High: 3 hrs. 45 min.; Low: 6 hrs.

BATTERIES INCLUDED: 6



*Dragonfly PRO features Cerakote, a ceramic polymer coating that offers industry-leading durability, hardness, scratch resistance, corrosion resistance, and chemical resistance.





Our RDH Model features the Orascope line of RDH Elite™ telescopes, which were designed in collaboration with dental hygienists, specifically for dental hygienists.

Our entry-level NEO model is available in a single color and supports 4 popular magnification powers. This is the preferred loupe for dental students and medical residents.



MAGNIFICATION

Traditional: 2.5x
Ergo: 3.0x

LIGHT INTENSITY 78 Lumens

BATTERY RUN TIME (Set of 2) 3 hrs. 45 min.

BATTERIES INCLUDED: 4



CHAMPAGNE



LILAC



MIDNIGHT

MAGNIFICATION

Traditional: 2.5x, 3.0x, 3.5x
Ergo: 3.5x

LIGHT INTENSITY 78 Lumens

BATTERY RUN TIME (Set of 2) 3 hrs. 45 min.

BATTERIES INCLUDED: 4



MIDNIGHT

PERFORMANCE
UNBOUND

Schedule your free in-office demo at Orascope.com/dragonfly



DRAGONFLY™
ORASOPTIC

PERFORMANCE UNBOUND



ORASOPTIC™

| **SUPERIOR VISUALIZATION**

6650 Ronald Reagan Avenue | Madison, WI 53704 | 800.369.3698 | [Orasoptic.com](https://www.Orasoptic.com)

United Kingdom | +44 1733 315203 | CustomerCareUK@Orasoptic.com

Nordics | +46 (0) 2 008 83009 | customercare.europe@orasoptic.com